



*Solutions For Life*



**Powered by...**



**Quantumrehab.co.uk**

# Stealth Cushion Range



## SIMPLICITY G



### Features:

- Contoured, high-density molded foam base
- Coccyx and sacral relief
- Moisture resistant and easy to clean Stealth-Tek incontinence cover.
- Sizes available - 16" - 24"W x 16" - 20"D
- Mesh cover optional



## SPECTRUM FOAM SPP



### Features:

- Lateral and pre-ischial shelf contoured, high-density molded foam base
- Multi-layer, viscoelastic foam insert under at risk bony areas
- CoolCore® technology cover standard.
- Urethane liner for incontinence, 4-Way Stretch Poly liner optional
- Sizes available - 16" - 32"W x 16" - 28"D
- Custom sizes available
- Optional Stealth-Tek or Mesh cover.



## SOLUTION SPP



### Features:

- Lateral and pre-ischial shelf contoured, high-density molded foam base
- Twin-cell, aqueous gel insert under at risk bony areas
- CoolCore® technology cover standard.
- Urethane liner for incontinence
- Sizes available - 10" - 24"W x 10" - 24"D
- Custom sizes available
- Optional Stealth-Tek or Mesh cover.



## SPECTRUM GEL SPP



### Features:

- Lateral and pre-ischial shelf contoured, high-density molded foam base
- Twin-cell, lightweight fluid insert under at risk bony areas
- CoolCore® technology cover standard.
- Urethane liner for incontinence
- 4-Way Stretch Poly liner optional
- Sizes available - 14" - 24"W x 16" - 24"D
- Custom sizes available
- Optional Stealth-Tek or Mesh cover.



### TRU-COMFORT 2 SPP



Convenient side pouch on both sides of the cushion with easy pull rubber loop zipper handles. Pouch is large enough to fit a larger format mobile phone.

### Features:


- Deep contoured, high-density molded foam base
- CoolCore Technology cover standard.
- Easy pull rubber loop zipper handles
- Convenient side pouches on both sides of the cushion
- Sizes available - 16" - 22"W x 14" - 24"D
- Custom sizes available
- 4-Way Stretch Poly liner optional

## Powered by...



### What is CoolCore®

- Chemical Free material
- UPF/Flammability rated
- 4-Way Stretch nylon interlock with Elastane
- Patented combination of chemical-free blends that deliver 3 distinct functions: **Moisture Wicking to Moisture Transportation to Moisture Evaporation.**
- **CoolCore® Technology is exclusive to Stealth Products.**



### Year Round Temperature Control

CoolCore fabrics wick moisture away from the body in every direction keeping you dry and comfortable no matter what season it is.

### Breakthrough Technology

Research driven, patented fabrics from Coolcore are chemical-free and yield remarkable benefits. They cool better, last much longer and offer protection from the sun's harmful rays (UPF).

## Stealth Cushion Specifications

Cushion Name	Weight Capacity	Weight of Cushion	Cushion Height	Width Range	Depth Range
<b>Simplicity G</b>	25 Stone	1 k.g	2" (5.08 cm)	16" - 20" (40.64 cm - 60.96 cm)	16" - 20" (40.64 cm - 60.96 cm)
<b>Solution SPP</b>	25 Stone	2.44 k.g	3" (7.62 cm)	10" - 24" (25.4 cm - 60.96 cm)	10" - 24" (25.4 cm - 60.96 cm)
<b>Spectrum Foam SPP</b>	46 Stone	1.36 k.g	3" (7.62 cm)	16" - 28" (40.64 cm - 81.28 cm)	16" - 28" (40.64 cm - 81.28 cm)
<b>Spectrum Gel SPP</b>	25 Stone	1.58 k.g	3" (7.62 cm)	16" - 24" (35.56 cm - 60.96 cm)	16" - 24" (35.56 cm - 60.96 cm)
<b>Tru Comfort 2 SPP</b>	32 Stone	1.81 k.g	4.5" (11.43 cm)	14" - 24" (40.64 cm - 55.88 cm)	14" - 24" (40.64 cm - 55.88 cm)



*Solutions For Life*

**Quantumrehab.co.uk**